



# Best<sup>®</sup> medical *international*

## Brachytherapy Needles



**Best Medical can  
provide any standard  
or customized  
configuration made  
to your exact  
specifications, sterile  
or non-sterile,  
shipped within 24  
hours — worldwide!**

## **Contact us:**

**Best Medical International  
7643 Fullerton Road  
Springfield, Virginia 22153 USA**

**Tel: 703-451-2378 • 800-336-4970**

**Fax: 703-451-5228**

**[www.bestmedical.com](http://www.bestmedical.com) • [www.teambest.com](http://www.teambest.com)**



## Mission Statement

*Everyone deserves the best health care. Our goal is to work with medical professionals to provide the Best products and services. As a leader in innovative technology, our mission is to uphold our reputation for excellence in the health care field by developing, manufacturing, and delivering cost-effective, high quality products to benefit patients around the world.*

## Corporate Values

*Fulfilling our mission requires a dedicated employee team—one that is committed to uncompromised quality and safety, and fostering a sense of pride in the company.*

## Best employees will:

- *Honor the dignity and uniqueness of each person, including fellow employees, customers, distributors, and vendors, through their actions, words, and behavior. They will treat others with respect, as they would wish to be treated.*
- *Embrace the spirit of cooperation and teamwork for which Best Medical is known, an environment which supports and recognizes effort and results.*
- *Promote and practice open and honest communication at all times.*
- *Listen to others and remain open to suggestions at all times, as we focus on achieving excellence in our work.*
- *Attract and retain leaders with exceptional abilities and integrity.*
- *Ensure all resources are managed in a responsible and effective manner.*

# Best ADVANTAGE

- Best Medical International is the **ONLY** supplier of Iridium-192 and Gold-198 seeds.
- We are the first and only provider of Iridium-192 in nylon ribbons in a wide range of activities.
- We are the first and only company to provide custom service for Gold-198 permanent implants, as well as a wide range of activities.
- We are the first and only company that will call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

---

# Catalog 2013 **Table of Contents**

Best Mick Type Needles ..... 1

Best Prostate Seeding Needles ..... 2

Prostate Stabilization Set Needles ..... 3

Needles ..... 4



*Best*<sup>®</sup> *medical international*


---

White Hub  
17 Gauge

Yellow Hub  
18 Gauge

The Mick type needle's obturator is provided with a sharp trocar point while the echogenic canula tip is ground-beveled to facilitate the penetration of the skin.

The Mick type needle's hub canula is "honed" for smooth transfer of seeds from the applicator into and through the needle.



### Best Mick Type Needles

Part #	
10203	17 Gauge White Hub
10209	18 Gauge Yellow Hub

### Best Mick Type Needles

- Sterile/disposable
- Mirror finish stylet
- Siliconized/calibrated
- Flared canula hub
- Improved sharpness
- Color coded hubs
- Trocar point
- Available 17 and 18 gauge
- Rigid stylet
- Echogenic bevel

### Packaged, Sterile and Disposable

- 25 Needles per carton
- 4 Pouches with 5 needles
- 5 Pouches with 1 needle

### Specifications (For Single Use Only)



## Best Prostate Seeding Needles

Best Prostate Seeding Needles	
Part #	
10215	Best C.P. Type Round Hub Needles
10201	Best C.P. Type Square Hub Needles
10197	Best Square Hub Needles
10199	Best Round Hub Needles

## Best Prostate Seeding Needles

- Sterile/disposable
- Low price
- Echogenic bevel
- Siliconized
- Calibrated stylet
- Improved sharpness
- Bevel point
- Available 17 and 18 Gauge

## Packaged, Sterile and Disposable


- 25 Needles per carton
- 5 Pouches with 1 needle
- 4 Pouches with 5 Needles each

## Specifications (For Single Use Only)

“A” and “B” bevel needles are available.







## Prostate Stabilization Set Needles

### Part #

10207

18 Gauge x 20 cm  
Sterile; Box of 10 Sets;  
2 Needles per Set

Improved Seed Distribution by Means of Prostate Gland Immobilization

### Prostate Stabilization Set

The Stabilizing Needle (patented) is designed to prevent prostate movement during the insertion of prostate seeding needles. Immobilization is accomplished by deployment of a “side barb” in the prostate tissue and then securing a “lock spring” against the proximal side of the template. Two stabilizing needles must be used to adequately stabilize the prostate gland.

### Features

- Two-part hub design
- Locking side barb
- Sharp bevel-point needle tip
- Echogenic tip
- 10 Pouches with 2 needles

### Packaged:

18 Gauge x 20 cm sterile; box of 10 sets; 2 needles per set

Needles	
Part #	
10171	Brachytherapy Needles Round Hub 18 Gauge x 10 cm with Wax in a Pouch
10180	Brachytherapy Needles Square Hub 18 Gauge x 20 cm with Wax in a Pouch
10183	Brachytherapy Needles Square Hub 17 Gauge x 20 cm without Wax in a Pouch
10222	Brachytherapy Needles Square Hub 17 Gauge x 20 cm with Wax in a Pouch
10191	Brachytherapy Needles Round Hub 18 Gauge x 20 cm with Wax in a Pouch
10195	Best Applicator Needles with "B" Bevel Blunt Stylet in a Pouch
10207	Prostate Stabilizing Needles 18 Gauge x 20 cm in a Pouch
10211	Brachytherapy Needles Round Hub 18 Gauge x 20 cm Pink Stylet with Wax in a Pouch
10213	Brachytherapy Needles Square Hub 18 Gauge x 20 cm Pink Stylet with Wax in a Pouch
10220	Brachytherapy Needles Round Hub 18 Gauge x 30 cm with Wax in a Pouch
10221	Brachytherapy Needles Round Hub 18 Gauge x 30 cm without Wax in a Pouch
10224	Needles Round Hub 16 Gauge x 20 cm without Wax (sterile)
10226	Needles Square Hub "B" Bevel 18 Gauge x 20 cm with Wax (sterile)
10229	Needles Round Hub "B" Bevel 18 Gauge x 20 cm with Wax (sterile)
10232	Brachytherapy Needle Metal Hub Zebra 18 Gauge x 20 cm
10211	Needles Square Hub 18 Gauge x 20 cm with Pink Stylet without Wax
5006-0033	18 Gauge x 20 cm Long Best Custom Stylet Applicator Needle (dittoli type)



**Best Medical International, Inc.** as author/distributor, believes the technical information and statements contained in this catalog to be reliable and correct at the time of publication. However, Best Medical International, Inc., does not assume any liability for statements, information or recommendations based on/or made in reliance on the information contained in this publication.

Data provided herein are for informational purposes only and are not intended for clinical use of any nature. Use of these data in such a manner is the sole responsibility of the user, and the author/distributor of the information makes no claim or inference as to its clinical suitability and/or accuracy. Best Medical International, Inc. assumes no warranty or liability for use of this information or products. Designs and specifications are subject to modification and/or change(s).

© 2013 Best Medical International, Inc.

## Terms and Conditions

**Ordering:** Orders may be placed with a purchase order number and transmitted by FAX, phone or e-mail.


**Payment:** Within 30 days.

**Freight:** FOB origin, freight prepaid and charged back through invoice.

**Returns:** Unused accessories may be returned to Best Medical International, Inc. within 15 days if found unsuitable for a particular application. No return or authorization number is necessary. A note stating the reason for return, name, and address of hospital will expedite processing of any credit.

**Note:** Best Medical International, Inc. assumes no liability for consequential damages of any kind for the material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer.

Best Medical International, Inc., reserves the right to modify the dimensions of accessories/components as well as the prices.



# Are you Searching for the Greatest Healthcare Company in the World?

## Look to TeamBest®!

World population is increasing by more than 100 million per year. To meet the needs of our aging and growing population, TeamBest® companies are committed to making quality healthcare affordable and accessible globally.

**Ready to Serve Globally for All Your Needs**

- Arplay Medical
- Best Cyclotron Systems
- Best Medical Africa
- Best Medical Asia
- Best Medical Canada
- Best Medical China
- Best Medical Engineering
- Best Medical Europe
- Best Medical International
- Best Medical Latin America
- Best Medical Middle East
- Best Medical Parts
- Best Medical Research
- Best Medical Russia
- Best NDT
- Best NOMOS
- Best Particle Therapy
- Best Theratronics
- Best Vascular (Novoste)
- CNMC Company
- Genetherapy Pharmaceuticals
- Huestis Diagnostics
- Huestis Medical and more...

**For more information, visit [www.teambest.com](http://www.teambest.com)**

© 2013 Best Medical International, Inc.



Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
703 451 BEST (2378) 800 336 4970 [www.bestmedical.com](http://www.bestmedical.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA