

## T-3 Mineral Filled Thermometer

The Model T-3 is a pocket-size, wide range, total-immersion, mineral-filled glass thermometer. The red disazo-type non-toxic dye is easily viewed against the white background. It is supplied in an aluminum protective casing with a ring on top and a pocket clip.

Range..... -10° C to 110° C  
 Accuracy..... ±0.25° C  
 Resolution ..... 1° C/div. 0.62 mm/° C  
 Length..... 14.9 cm (5.5 in)  
 Weight..... 23 g (0.8 oz)

## T-5 PerformaTherm™ Filled Thermometer

The Model T-5 is a precision, total immersion, proprietary non-toxic biodegradable liquid filled thermometer designed specifically to cover limited range of temperatures encountered in a typical radiation therapy department, thereby offering excellent resolution, combined with a short length. Its blue liquid column is easily visible against the white background. Includes a NIST traceable calibration certificate.

This thermometer may be ordered with an open-face, brass, armored case as Model T-5A.

Range..... -1° C to 27° C  
 Accuracy..... ±0.2° C  
 Resolution ..... 0.2° C/div., 5 mm/° C  
 Length..... 22.9 cm (9 in)  
 Weight..... 34 g (1.2 oz)



**Model T-3**



**Model T-5**

